

Miniature Raman Spectrometer

ATR2500

Features

- Full free space design, super high sensitivity;
- Ultra-high sensitivity detector;
- Ultra-low noise circuit;
- Ultra-light, ultra-small, ultra reliable;
- Powerful spectral analysis software;
- Eliminate fluorescent background;
- Peak search and display;
- USB 2.0;
- Friendly man-machine interface.

Application

- Research of Research Institute
- Nanoparticles and new materials
- Bioscience
- Forensic Medicine
- Material science
- Medical immunoassay
- Agriculture and food identification
- Gems and inorganic minerals ID
- Environmental science

Description

ATR2500 is a Raman spectrometer developed by Optosky for more than 20 years. After 5 years, it has developed a brand-new, optimized and designed high-sensitivity Raman spectrometer with breakthrough characteristics. It has ultra-small and ultra-light, High resolution, high sensitivity, high reliability, etc. ATR2500 adopts Optosky's latest full free space optical path technology, which increases the Raman signal collection efficiency by nearly 4 times, thereby increasing the sensitivity by 4 times.

The ATR2500 Raman spectrometer is very suitable for laboratory scientific research. It is small in size, high reliability, easy to measure, and the detection results are accurate and reliable. The excellent low stray light design of ATR2500 makes it easy to use. The multi-function software randomly distributed by ATR2500 has been strictly tested by hundreds of scientists around the world and collected their improvement opinions. After nearly a hundred versions of updates, the function is very complete and stable, which is very suitable for the development of Raman research.



1. Selection Guide

Mode	Spectral Range (cm ⁻¹)	Resolution (cm ⁻¹)	Excitation Wavelength (cm ⁻¹)
ATR2500-R27	250-2700	3~6	785
ATR2500-R40	150-4000	7~10	785
ATR2500-532	100-4100	8~12	532
ATR2500-1064	200-2600	10~15	1064

2. parameter

electronic Parameters	
Interface	USB 2.0
Integration time	4ms - 120s
Voltage	DC 5V±5%
Work temperature	-10~45 °C
Work humidity	< 95%
Dimension(L*W*H)	Without probe: 119.2×89×35 mm With probe: 139×89×35 mm
Weight	390 g
Reliability	
Spectral reliability	$\sigma/\mu < 0.5\%$ (COT 8 hours)
Temperature reliability	Spectral shift ≤ 1 cm ⁻¹ (10-40 °C)
Spectral intensity	$< \pm 5\%$ (in 5 ~ 40 °C)
SNR	>1300:1 (918 cm ⁻¹ of Acetonitrile , 4sIntergation, 130mW)
Detector	
Model	Ultra-sensitive linear array detector
Spectral range	200-1100 nm
Effective pixels	2048 pixels
Dynamic range	50000: 1
Laser	
Center wavelength	785±0.5nm
Half width	0.08 nm
Maximum output power	≥ 300 mW , the actual output power software can be set

Minimum power output adjustment	1mW
Power stability	$\sigma/\mu < \pm 0.2\%$
Raman probe	
Rayleigh Scatter Resistance	6 mm
Operating Distance(OD)	OD>8
NA	0.3
Aperture	7mm

3. ATR2500 Spectrum




Figure 1 Raman spectrum curve test results; sample: acetonitrile, laser power: 130mW, measurement integration time: 4000ms.




Figure 2 Spectral intensity repeatability test, the result is 0.06%, the spectral intensity stability is good.




Figure 3 Temperature drift test, -10~40°C wavenumber drift

4. ATR2500 Dimensions




Figure 4 ATR2500 Dimensions

5. ATR2500 Series

Optical Parameters	ATR2500	ATR2500-532	ATR2500-1064
Excitation wavelength	785±0.5nm	532±0.5nm	1064±0.5nm
Maximum laser power	130	70	500
Laser linewidth	0.08	0.1	0.08
Spectral range	150-4000 or 250-2700	100-4100	200-2600
Resolution	7-10 or 3/4-6	8-12	10-15