

# INOVATION CARE AND HEALTH



Wehere Medical Technology Co.,Ltd  
*High-quality and High-output*



# INTRODUATION

Hangzhou Wehere Medical Technology is a company on the medical health products manufacturing and trading in China.

As a professional medical supplies over twenty years, our products have exported globally to Europe, Asia, America, etc., and help customers to succeed in their markets.

We have three plants to produce different kinds of medical products, and some of them can be customized.

On the basis of production management MDR and ISO13485 in our owner factories, guided by brand strategy, under the direct participation of American companies, we are committed to becoming an important member of the global medical and healthcare industry.

We are here to support your medical and healthcare needs.



# Urine Bag Animal Drainage Collection Bag

WH-0001(300ML)



Drain valve: push-pull valve

Configuration: includes two fixing straps and locating buckles

Backflow prevention

The connecting device is provided with a luer connector

WH-0002(350ML)



Drain valve: Tvalve

Configuration: includes two fixing straps and locating buckles

Backflow prevention

The connecting device is provided with a luer connector

WH-0003(300ML)



Drain valve: push-pull valve

Configuration: includes two fixing straps and locating buckles

Backflow prevention

The connecting device is provided with a luer connector

WH-0004(150ML)



Drain valve: push-pull valve

Configuration: includes two fixing straps and locating buckles

Backflow prevention

The connecting device is provided with a luer connector

# Accessories



Model: Fr6, Fr8, Fr10  
Configuration: Fixed Balloon and Guide Wire

Model: 1.0\*130MM/1.3\*130MM  
Configuration: PU tube and guide wire



WH-0005  
2-way Silicone Catheter



Material: PVC

Model: 4.0\*150CM/4.0\*200CM  
Configuration: Luer Connector and Cap



WH-0006  
Hydrophilic Ureteral Access sheath



WH-0007  
Medical Urine Bag Connector

WH-0008  
PU Spiral Tubes



# Different Types of Pet Urine Pad



WH-0009

Pet Pee Pad (Ultrasonic Series)

## Product composition:

Hydrophilic non-woven, Dust-free paper, Japan SAP polymer, PE coated non-woven fabric

## Product size:

S: 33\*45 M: 45\*60 L:60\*60 XL:60\*90

## Super Absorption & Altrasonic Seamless Bonding Technology

This Technology can enhance water resistance without sewing threads and pinholes, so as to meet the requirements of watertight, airtight, anti-virus and anti-bacteria.

## Product composition:

Hydrophilic soft non-woven, Sanitary paper, Imported fluff pulp, SAP polymer, PE leakproof bottom film

## Product size:

S: 33\*45 M: 45\*60 L:60\*60 XL:60\*90

**Maternal and Infant Standards Hydrophilic non-woven Fabric, Soft as Clouds**



WH-0010

Pet Pee Pad (Cloud Series)

# Different Types of Pet Urine Pad



WH-0011  
Pet pee Pad (Silk Series)

## Product composition:

Hydrophilic non-woven, Sanitary paper, Imported fluff pulp, SAP polymer, PE leakproof bottom film

## Product size:

S: 33\*45 M: 45\*60 L:60\*60 XL:60\*90

## Thin and Absorbent & Smooth to the Touch

## Product composition:

Non-woven, Sanitary paper, Bamboo carbon fiber, Imported raw wood pulp, SAP polymer, PE

## Product size:

S: 33\*45 M: 45\*60 L:60\*60 XL:60\*90

## Beodorization and Odor Absorption



WH-0012  
Pet pee Pad (Panda Series-  
Bamboo Charcoal Fiber)

# Different Types of Pet Urine Pad



WH-0013  
DIY Pet Pee Pad (Flora Series)

## Product composition:

Hydrophilic soft non-woven, Sanitary paper, Deodorant factor, Imported wood pulp, SAP polymer, PE leakproof bottom film

## Product size:

S: 33\*45 M: 45\*60 L:60\*60 XL:60\*90

## Thin Full of Deodorant Factor & Pleasant Smell


## Product composition:

Hydrophilic non-woven, Sanitary paper, Imported wood pulp, SAP polymer, Super strong anti-leakage bottom film

## Product size:

S: 33\*45 M: 45\*60 L:60\*60 XL:60\*90

**Used as Cushion Stretcher, the Maximum Weight of 170KG**



WH-0014  
Pet Pee Pad (Apollo Series)

# Product Parameter



## Product Name:

Fast-absorbent Pet Pad (Sponge series)

## Product Composition:

PET non-woven Japan SAP polymer

## Product Size:

S: 25\*30 M: 30\*40 L: 40\*50



WH-0015  
Thickening Type



WH-0016  
Activated Carbon Type



WH-0017  
Normal Type



WH-0018  
Degradable Type



# Solve Your Pet Problems

Soft, Comfortable, Long-lasting, Dry, Waterproof, Urine repellent,  
Instant absorption, Leak proof

Hydrophilic soft non-woven fabric

Maternal and infant grand sanitary paper

Polymer imported from Japan

Imported fluff pulp

PE anti-leakage base film

Unique embossed flow guide layer design

Lock up the urine quickly

Reverse osmosis rejection



## Why are we so dry?

Six-layer three-dimensional water locking process

1. Coated sanitary paper
2. Soft breathable layer is comfortable and skin-friendly
3. Imported fluff pulp instant urine absorption and non-diffusion
4. Polymer SAP absorbs several times the volume of water
5. Maternal and infant grade toilet paper wrapped cotton core
6. PE composite substrate prevents urine leakage

# Handheld Pain Relief Device



## WH-0019

1. Pet animal pain relief, Knee arthritis, Rheumatoid Arthritis, Athletic system injuries, Soft tissue injuries, Relieve muscle and joint pain, Surgical incisions, Tendon and ligament injuries, Traumatic injuries.

2. Diminish inflammations, Reduce swelling, Speeding up wound healing, Ulcers, Osteoarthritis, Acupuncture, Rehabilitation therapy

The veterinary industry has deemed pet laser therapy safe and effective. It can effectively treat diseases, injuries and conditions such as tissue injuries (including strains and sprains) and arthritis. This non-invasive, drug-free treatment is used to treat inflammatory conditions.


More recently, it has been used in conjunction with traditional medical treatments to treat soft tissue or tendon injuries and arthritis. It can also be used to encourage wound healing.



Before



After



We are excited to offer laser therapy as a non-invasive, drug free treatment for healing.

# Operation Instructions:

1.Preparation before treatment: Use warm water to clean skin first; please keep away from hair and black nevus of your skin; prevent laser energy burning skin.

2.Turning on: long pressing the power key for 2 seconds, the screen displaying.

3.Selecting treatment time: pressing "T" button, there are 10 minutes/ 15 minutes/ 20 minutes/ 25 minutes/ 30 minutes optional.

4.Selecting treatment power: pressing the "P" button, there are four intensity can be adjusted.

5.Starting treatment: Again short press the "ON/OFF" to start treatment, you will see the laser probes are emitting the laser light. Put the device to the pain area.

6.Turning off: long pressing On/Off button, indicated light turning out, the screen no display.

7.Charging: the instrument with lithium battery inside , when instrument have the invoice of "di,di,di" ,please recharge,the each recharging time will be 8~10 hours : when instrument is charging ,please do not use.

## Notice:

- 1.Clear the skin before treatment;
- 2.Avoid the sweat, anything in black color like hairs and mole;
- 3.Any uncomfortable, please just remove the device to stop treatment.



## TREATMENT PRINCIPLE:

### 1. Pet animal arthritis


LLL T 1s a non - invasive light source treatment. It is also called photobiology or biostimulation. LLL T is believed to affect the function of cells , accelerate connective tissue repair and act as an anti inflammatory agent. Lasers with 650nm and 808nm are used in the treatment animal arthritis.

### 2. Pet Wound healing

Low level laser stimulate the cells to increase the production of certain enzymes which affect cell proliferation or cell division, thus increasing the healing speed of a wound

### 3. Animal Sports injures rehabilitation


LLL T use light photons to irradiate an area of injury to relief pain and promote animal natural heal-ing process. It reduce inflammation and encourage cell reproduction, offering pain relief and treat injures.





## Application:

1. Pet animal pain relief, Knee arthritis, Rheumatoid Arthritis, Athletic system injuries, Soft tissue injuries, Relieve muscle and joint pain.
2. Diminish inflammations, Reduce swelling, Wounds & Ulcers, Osteoarthritis, Acupuncture, Rehabilitation therapy.



Before



After

## Laser Therapy – A New Solution for Healing

A beagle dog is shown looking towards a handheld laser therapy device. The device is white with a black circular control panel featuring a digital display showing '88:88' and several buttons. Red lines representing the laser beam are shown emanating from the tip of the device. A speech bubble from the dog contains text about the benefits of laser therapy.

*We are excited to offer laser therapy as a non-invasive, drug free treatment for healing.*

**Laser Therapy has been successful in treating post surgical pain and many acute or chronic conditions. Laser Therapy may be beneficial if your pet is in pain, has inflammation or has a wound.**


# Veterinary Low Level Laser



## WH-0020 Therapy Pad Wrap

### ADVANTAGES

- ◆ Enhance circulation
- ◆ Speed the healing process
- ◆ The release of endorphins
- ◆ Vasodilation can increase blood flow ultimately increasing oxygen and speeding up the healing



Allow the muscles to relax

Reduce pain and swelling

Decrease nerve sensitivity

In addition, laser therapy does not have any negative side effects and no sedation is required. You also don't need to worry about your pet missing any fur as cold laser therapy doesn't require the removal of fur in order to properly complete treatment.

Veterinary portable pain relief pad is a low level laser therapy (LLL T) device, also known as a soft laser. LLL T is a light source treatment that generates a coherent light of a single wavelength.

It emits no heat, sounds or vibration. The laser's energy is used to penetrate the skin and stimulate the body's natural healing processes. By this way, every cell in the treatment area received the energy activate deep tissue repair, Laser irradiation can diminish inflammation to edema fluid, accelerate the secretion of analgesic sub-stance and reduce the sensitivity of the peripheral nerve

### HOW TO USE IT

Tips: Use this pad stick on below parts for different use. We suggest to use the device for 2 times per day, 15-30 minutes per time.

Do not continuously use the device more than 30 minutes per time to avoid the device damage.



# Hyperbaric Chamber Veterinary Type

WH-0021

| Model          |                | 1                                                                                | 2                               |
|----------------|----------------|----------------------------------------------------------------------------------|---------------------------------|
| Chamber Cabin  | Pressure       | 1.3/1.4ATA(4/5.8psi)                                                             | 1.5ATA(7psi)                    |
|                | Material       | Aluminum+Acrylic                                                                 | Stainless Steel +PC             |
|                | Cabin diameter | 60cm                                                                             | 70cm                            |
|                | Size(D*L)      | 125x73x86cm (49x29x34 inch)                                                      | 142x71 x93cm<br>(57x28x37 inch) |
|                | Weight         | 65kg                                                                             | 100kg                           |
| System Machine | Size/Weight    | 77x46x72cm (30x18x28 inch), 70kg                                                 |                                 |
|                | Flow of Air    | 72 Liter/min                                                                     |                                 |
|                | Flow of Oxygen | 5Liter/min /10 liter/min                                                         |                                 |
|                | Air Filter     | Dual                                                                             |                                 |
|                | Voltage        | 11 OVAC/220VAC $\pm$ 10%                                                         |                                 |
|                | Description    | All-in-one machine, a combine of air compressor, oxygen concentrator, and cooler |                                 |



## Application Fields:

Gastric dilation and volvulus, and the associated injuries.

Fracture repair, especially open fractures or those involving malunion or nonunion

Gastrointestinal resection and anastomosis, especially post-foreign body obstruction.

Amputations

Skin grafts












Hemilaminectomy

Ear canal ablation






# Rapid Test




## Easy to Use Easy to Read


- ◆ Take out swab
- ◆ Take samples
- ◆ Mixed
- ◆ Take mixed buffer
- ◆ Add 3-4 drops of buffer
- ◆ 5-10 mins
- ◆ Read results



Negative



Positive



Invalid

# Product List-Canine Rapid Test



| Product Name                            | Abbr.       | Specimen     | Format   | Specification |
|-----------------------------------------|-------------|--------------|----------|---------------|
| Canine Distemper Virus Antigen Test     | CDV Ag      | Secretions   | Cassette | 20T           |
| Canine Distemper Virus Antibody Test    | CDV Ab      | Serum/Plasma | Cassette | 20T           |
| Canine Parvo Virus Antigen Test         | CPV Ag      | Feces        | Cassette | 20T           |
| Canine Parvo Virus Antibody Test        | CPV Ab      | Serum/Plasma | Cassette | 20T           |
| Canine Influenza Virus Ag Rapid Test    | CIV Ag      | Secretions   | Cassette | 20T           |
| Canine Coronavirus Antigen Test         | CCV Ag      | Feces        | Cassette | 20T           |
| Canine Parainfluenza Virus Antigen Test | CPIV Ag     | Secretions   | Cassette | 20T           |
| Canine Adenovirus II Antigen Test       | CAV-II Ag   | Secretions   | Cassette | 20T           |
| Canine Adenovirus I Antigen Test        | CAV-I Ag    | Secretions   | Cassette | 20T           |
| Canine CRP Test                         | C-CRP       | Serum/Plasma | Cassette | 20T           |
| Canine Toxoplasma Antibody Test         | TOXO Ab     | Serum/Plasma | Cassette | 20T           |
| Canine Heartworm Antigen Test           | CHW Ag      | Serum/Plasma | Cassette | 20T           |
| Leishmania Canis Antibody Test          | LSH Ab      | Serum/Plasma | Cassette | 20T           |
| Canine Brucella Antibody Test           | C.Bru Ab    | Serum/Plasma | Cassette | 20T           |
| Ehrlichia Canis Antibody Test           | EHR Ab      | Serum/Plasma | Cassette | 20T           |
| Canine Leptospirosis Antibody Test      | Lepto Ab    | Serum/Plasma | Cassette | 20T           |
| Babesia Gibsoni Antibody Test           | BG Ab       | Serum/Plasma | Cassette | 20T           |
| Rabies Antigen Test                     | Rabies Ag   | Secretions   | Cassette | 20T           |
| Rabies Antibody Test                    | Rabies Ab   | Serum/Plasma | Cassette | 20T           |
| Lyme Disease Antibody Test              | Lyme Ab     | Serum/Plasma | Cassette | 20T           |
| Pregnancy Relaxin Test                  | RLN         | Serum/Plasma | Cassette | 20T           |
| Canine Giardia Antigen Test             | C-GIA Ag    | Feces        | Cassette | 20T           |
| CDV/CPIV Ag Combo Test                  | CDV/CPIV Ag | Secretions   | Cassette | 20T           |
| Canine Parvo/Corona Ag Combo Test       | CPV/CCV Ag  | Feces        | Cassette | 20T           |
| Canine Anaplasma Antibody Test          | C.ANA Ab    | Serum/Plasma | Cassette | 20T           |
| Canine Rotavirus Antigen Test           | ROTA        | Secretions   | Cassette | 20T           |
| CPV/CDV Antibody Combo Test             | CPV/CDV Ab  | Serum/Plasma | Cassette | 20T           |

# Product List-Canine Rapid Test



| Product Name                                                           | Abbr.              | Specimen                 | Format   | Specification |
|------------------------------------------------------------------------|--------------------|--------------------------|----------|---------------|
| Canine Distemper/Adeno Ag Combo Test                                   | CDV/CAV Ag         | Secretions               | Cassette | 20T           |
| Canine Parvo-Corona-Rota Virus Antigen Combo Test                      | CPV/CCV Rota Ag    | Feces                    | Cassette | 20T           |
| CPV/CCV/Giardia Combo Test                                             | CPV/CCV/Giardia Ag | Feces                    | Cassette | 20T           |
| Canine Distemper/Adeno Influenza Combo Test                            | CDV/CAV/CIV        | Secretions               | Cassette | 20T           |
| Canine Infectious Hepatitis/Parvo Virus/Distemper Virus IgG Combo Test | ICH/CPV/CDV        | Whole Blood/Serum/Plasma | Cassette | 20T           |
| Canine Ehrlichia/Anaplasma Combo Test                                  | EHR/ANAAb          | Whole Blood/Serum/Plasma | Cassette | 20T           |
| Ehrlichia/Lyme/Anaplasma Combo Test                                    | EHR/LYM/ANA Ab     | Serum/Plasma             | Cassette | 20T           |
| Ehrlichia/Lyme/Anaplasma/hearworm Combo Test                           | EHR/LYM/ANA/CHW    | Serum/Plasma             | Cassette | 20T           |
| Ehrlichia/Babesia/Anaplasma Combo Test                                 | EHR/BAB/ANA        | Serum/Plasma             | Cassette | 20T           |
| Ehrlichia/Babesia/Anaplasma/ hearworm Combo Test                       | EHR/BAB/ANA/CHW    | Serum/Plasma             | Cassette | 20T           |



# Product List-Feline Rapid Test



| Product Name                                                              | Abbr.          | Specimen                 | Format   | Specification |
|---------------------------------------------------------------------------|----------------|--------------------------|----------|---------------|
| Feline Panleucopenia Antigen Test                                         | FPV Ag         | Feces                    | Cassette | 20T           |
| Feline Infectious Peritonitis Antibody Test                               | FIP Ab         | Serum/Plasma             | Cassette | 20T           |
| Feline Infectious Peritonitis Antigen Test                                | FIP Ag         | Serum/Plasma/Feces       | Cassette | 20T           |
| Feline Coronavirus Antigen Test                                           | FCoV Ag        | Feces                    | Cassette | 20T           |
| Feline Leukemia Virus Antigen Test                                        | FeLV Ag        | Serum/Plasma             | Cassette | 20T           |
| Feline Immuno Deficiency Virus Antibody Test                              | FIV Ab         | Serum/Plasma             | Cassette | 20T           |
| Feline Giardia Antigen Test                                               | F-GIA Ag       | Feces                    | Cassette | 20T           |
| Feline Anaplasma Antibody Test                                            | ANA Ab         | Serum/Plasma             | Cassette | 20T           |
| Feline Toxoplasma Antibody Test                                           | TOXO Ab        | Serum/Plasma             | Cassette | 20T           |
| Feline Viral Rhinotracheitis Antigen Test                                 | FHV Ag         | Secretions               | Cassette | 20T           |
| Feline Calcivirus Antigen Test                                            | FCA Ag         | Secretions               | Cassette | 20T           |
| Feline Heartworm Antigen Test                                             | FHW Ag         | Serum/Plasma             | Cassette | 20T           |
| Feline panleucopenia Antibody Test                                        | FPV Ab         | Serum/Plasma             | Cassette | 20T           |
| Feline Coronavirus Antibody Test                                          | FCoV Ab        | Serum/Plasma             | Cassette | 20T           |
| Feline Herpes Virus Test<br>(Feline viral Rhinotracheitis Antigen Test)   | FHV Ag         | Secretions               | Cassette | 20T           |
| FIV Ab/FeLV Ag Combo Test                                                 | FIV Ab/FeLV Ag | Serum/Plasma             | Cassette | 20T           |
| Feline Herpes/Feline Calicivirus Virus Comb Test                          | FHV/FCA Ag     | Secretions               | Cassette | 20T           |
| Feline Panleukopenia/Herpes Virus/Calici Virus<br>IgG Antibody Combo Test | FPV/FHV/FCV    | Whole Blood/Serum/Plasma | Cassette | 20T           |



# Product List-Livestock Rapid Test



| Product Name                                                  | Abbr.       | Specimen     | Format   | Specification |
|---------------------------------------------------------------|-------------|--------------|----------|---------------|
| Porcine Rotavirus Antigen Test                                | PRV Ag      | Feces        | Cassette | 20T           |
| Swine Infectious Gastroenteritis Antigen Test                 | TGE Ab      | Feces        | Cassette | 20T           |
| Porcine Epidemic Diarrhea Virus Anti-IgA Test                 | PED IgA     | Serum/Plasma | Cassette | 20T           |
| Porcine Circovirus Antibody Test                              | PCV Ab      | Serum/Plasma | Cassette | 20T           |
| Porcine Trichinella Spiralis Antibody Test                    | PTS Ab      | Serum/Plasma | Cassette | 20T           |
| Classical Swine Fever Virus Antibody Test                     | CSFV Ab     | Serum/Plasma | Cassette | 20T           |
| Porcine Pseudorabies-gE Antibody Test                         | PRV gE Ab   | Serum/Plasma | Cassette | 20T           |
| Porcine Pseudorabies-gB Antibody Test                         | PRV gB Ab   | Serum/Plasma | Cassette | 20T           |
| Porcine PRRS Antibody Test                                    | PRRSV Ab    | Serum/Plasma | Cassette | 20T           |
| Swine Foot and Mouth Diseases Virus Serotype-O Antibody Test  | C.FMDV-O Ab | Serum/Plasma | Cassette | 20T           |
| Swine Foot and Mouth Diseases Virus Serotype-A Antibody Test  | C.FMDV-AAb  | Serum/Plasma | Cassette | 20T           |
| Newcastle Disease Virus Antigen Test                          | NDV Ag      | Secretions   | Cassette | 20T           |
| Avian Influenza Virus Antigen Test                            | AIV Ag      | Secretions   | Cassette | 20T           |
| Avian Influenza Virus H5 Antigen Test                         | AIV H5 Ag   | Secretions   | Cassette | 20T           |
| Avian Influenza Virus H7 Antigen Test                         | AIV H7 Ag   | Secretions   | Cassette | 20T           |
| Avian Influenza Virus H9 Antigen Test                         | AIV H9 Ag   | Secretions   | Cassette | 20T           |
| Bovine Foot and Mouth Diseases Virus Serotype-O Antibody Test | B.FMDV-O Ab | Serum/Plasma | Cassette | 20T           |
| Bovine Foot and Mouth Diseases Virus Serotype-A Antibody Test | B.FMDV-A Ab | Serum/Plasma | Cassette | 20T           |
| Bovine Brucella Antibody Test                                 | B.Brucella  | Serum/Plasma | Cassette | 20T           |
| Sheep Brucella Antibody Test                                  | B.Brucella  | Serum/Plasma | Cassette | 20T           |

# Product List-Livestock Rapid Test



| Product Name                                                            | Abbr.      | Specimen     | Format   | Specification |
|-------------------------------------------------------------------------|------------|--------------|----------|---------------|
| Bovine Virus Diarrhea Antibody Test                                     | BVDV Ab    | Serum/Plasma | Cassette | 20T           |
| Bovine Infectious Rhinitis Antibody Test                                | IBR Ab     | Serum/Plasma | Cassette | 20T           |
| Clostridium Perfringens Antibody Test                                   | CLP Ab     | Serum/Plasma | Cassette | 20T           |
| Clostridium Spoilage Antibody Test                                      | CLS Ab     | Serum/Plasma | Cassette | 20T           |
| Peste Petitis Ruminantis Antibody Test                                  | PPR Ab     | Serum/Plasma | Cassette | 20T           |
| African Swine Fever Virus Antibody Test                                 | ASFV Ab    | Serum/Plasma | Cassette | 20T           |
| African Swine Fever Virus Antigen Test                                  | ASFV Ag    | Secretions   | Cassette | 20T           |
| Foot and Mouth Diseases Virus Non-structural Protein 3ABC Antibody Test | FMDV(3ABC) | Serum/Plasma | Cassette | 20T           |




# Automatic Hematology Analyzer (3-part differential)

**WH-0022 (3000 Vet)**

## Species of Animal

Dog, cat, pig, cow,  
sheep, goat, horse,  
monkey, mouse, rabbit,  
pot-bellied pig, camel,  
buffalo and 2 self-  
programmed animals.



## TECHNICAL PARAMETERS

WBC, NEU#, MXD#, LYM#, NEU%,  
MXD%, LYM%, RBC, HGB, MCV,  
MCH, MCHC, RDW-CV, RDW-SD,  
HCT, PLT, MPV, PDW, PCT, P-LCR,  
P-LCC



## Fast Test

Use double-channel test mode with 60 T/H

## Small Sample Volume Needed

6  $\mu$ L of venous blood (whole blood mode)  
20  $\mu$ L of capillary blood (pre-dilution mode)

## Easy and Convenient Operation

Large LCD touch screen  
High-definition color monitor

## No Reagent Residue

No residue of hematology diluent or  
lyse reagent, cost-effective



# Automatic Hematology Analyzer (5-part differential)

**WH-0023 (5000 Vet)**

## Species of Animal

Dog, cat, pig, cow, sheep, goat,  
horse, monkey, mouse, rabbit,  
pot-bellied pig, camel, buffalo  
and 2 self-programmed animals.



## Principles

Laser scatter in flow cytometry  
Classic electrical impedance method  
Mature photoelectric colorimetry

## More Accurate Test

Digital filtering, coincidence correction in blood  
cell counting  
Use optical method to detect Basophil  
Curve fitting, floating threshold technology  
The smallest gem hole in the industry with high  
sensitivity Intelligent scatter classification  
technology

## TECHNICAL PARAMETERS

WBC, NEU#, LYM#, MON#, EOS#, BAS#, NEU  
%, LYM%, MON%, EOS%, BAS%, RBC, HGB,  
HCT, MCV, MCH, MCHC, RDW-CV, RDW-SD,  
PLT, MPV, PDW, PCT, P-LCR, P-LCC





# Pet Bed List -Care and Warmth



WH-0024



WH-0025



WH-0026



WH-0027



WH-0028



WH-0029



WH-0030



WH-0031



WH-0032



WH-0033



WH-0034



WH-0035

# Pet Bed List -Care and Warmth



WH-0036



WH-0037



WH-0038



WH-0039



WH-0020



WH-0041



WH-0042



WH-0043



WH-0044



WH-0045



WH-0046



WH-0047

# Pet Bed List -Care and Warmth



WH-0048



WH-0049



WH-0050



WH-0051



WH-0052



WH-0053



WH-0054



WH-0055



WH-0056



WH-0057



WH-0058



WH-0059



# Orthopedic Power Tools System



Canulate Drill(B140)  
WH-0060




Multi-functional Power Drill Saw System(B100)  
WH-0061



Canulate Drill(B-D140)  
WH-0062



Multi-functional Power Drill Saw System(E-B100)  
WH-0063



TPLO Saw  
WH-0064



Oscillating Bone Saw(B150)  
WH-0065

# Portable X-ray Machine



Appearance C design, ABS engineering, plastic manufacturing, beautiful light and firm  
The overall size is 51 cm\*36cm\*2cm, which is convenient to carry  
The screen displays high-definition black and white images, without spots  
No darkrooms is needed for operating indoors or outdoors More optional functions  
One-key switch machine



## Working Principle/ Organization Structure

It uses X-ray to penetrate the measured object. The image of the internal structure of the object being measured is displayed on the LCD screen in front of the X-ray receiver.

It consists of high vol cage generator, X-ray tube, X-ray receiver, beam, frame, adapter etc.

## Scope of Application

It's mainly used in veterinary (smaller & thinner parts of animal) for X-ray diagnosis.

Scientific research and testing occasions

Other Non-destructive testing applications



## After-sales Service

After the customer purchases, we promise that the Portable X-ray Machine will be repaired freely for one year.

It can also be tailored to the needs of customers.



# Medical Equipment

Altrasound  
Machine WH-0066



Vet Electrocardiograph  
WH-0067



Veterinary Anesthesia  
Machine WH-0068



Veterinary Biochemical  
Instrument WH-0069



Therapy Cage  
WH-0070



Hospital Cage  
WH-0071



# Medical Equipment



Electric Beauty  
Table WH-0075



Operating Table  
WH-0072



Sink  
WH-0073



Disposal Table  
WH-0074



Automatic Pet Dryer  
Room WH-0076

# CONTACT US

#701, Block C, Wanda Plaza, Gongshu, Hangzhou,  
310015, China

Tel: +86-15381022017

Email: [della\\_liu@weheremed.com](mailto:della_liu@weheremed.com)

[www.weheremedical.com](http://www.weheremedical.com)

